**MANUAL DE INSTRUCCIONES
OPERATING MANUAL**

**THERMOVISC 100-F8
THERMOVISC 200-F8**

**ULTRATHERMOSTATOS DE CIRCULACIÓN *CIRCULATION*
ULTRATHERMOSTAT**

INFORMACIÓN GENERAL

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GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact with your supplier or FUNGILAB, S.A.'s technical service.*
- 5) **IMPORTANT! FUNGILAB, S.A. WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.*

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CONTENTS

GENERAL INFORMATION
CONTENTS
PACKING LISTS
TECHNICAL FEATURES
INSTALLATION
OPERATION
SPARE PARTS
GUARANTEE
COMMUNICATION RS-232
MAINTENANCE
"EC" CONFORMITY DECLARATION

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PACKING LIST

The standard equipment consist of the following components:

Pieza	Piece	Código / Code	Cant. / Qty.
Thermovisc 100-F8	<i>Thermovisc 100-F8</i>	B100F8	1
Thermovisc 200-F8	<i>Thermovisc 200-F8</i>	B200F8	1
Manual de instrucciones	<i>Instruction manual</i>		1

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TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Código / <i>Code</i>		B100F8				B200F8			
Capacidad litros / <i>Capacity litres</i>		8							
Rango de regulación / <i>Regulation range</i>		-10 / +100				-30 / +100			
		-10°C	0°C	+20°C	+100°C	-30°C	0°C	+20°C	+100°C
Estabilidad °C / <i>Stability °C</i>		±0,1	±0,1	±0,1	±0,1	±0,1	±0,1	±0,1	±0,1
Homogeneidad		±0,5	±0,1	±0,1	±0,1	±0,1	±0,1	±0,1	±0,1
Error de consigna		±0,25	±0,25	±0,25	±1	±0,5	±0,25	±0,25	±0,1
Resolución / <i>Resolution</i>		0,1	0,1	0,1	0,1	0,1	0,1	0,1	0,1
Medidas abertura de la cubeta cm. / <i>Tank opening dimensions cm.</i>	Ancho / <i>Length</i>	20							
	Fondo / <i>Width</i>	14							
Medidas exteriores cm. / <i>Overall dimensions cm.</i>	Alto / <i>Height</i>	26							
	Ancho / <i>Length</i>	66							
	Fondo / <i>Width</i>	44							
Bomba / <i>Pump</i>	Presión / <i>Pressure</i>	150 mbar							
	Caudal H ₂ O / <i>H₂O flow</i>	12 l/min.							
Consumo W. / <i>Consumption W.</i>		1150				1460			
Potencia compresor H.P. <i>Compressor power H.P.</i>		1/8				3/8			
Peso Kg. / <i>Weight Kg.</i>		28				30			

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
FILLING THE BATHS

For correct use of Immersion thermostats, first it is necessary that the tank to be used is filled with liquid up to adequate level, i.e. to at least cover the heating element and circulation pump.

Otherwise the thermostat does not work by float effect.

The liquid or mixture to be used must be adequated to the model of thermostat and temperature range.

As a guidance, here are some recommended liquids, however, other similar liquids can be used:

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INSTALLATION

Place the thermovisc on a flat, horizontal, level surface, trying to leave a space at the back of the apparatus of 10 cm.

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Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.

Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the apparatus.

Do not use the apparatus if it is NOT earthed.

If you have to change the plug bear in mind the following:


Blue cable: Neutral.
Brown cable: Phase.
Yellow/green cable: Earth.

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OPERATION OF THE PUMP WITHOUT EXTERNAL RECIRCULATION :

Never use the pump without the connection tube (see photo).

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OPERATION OF THE PUMP WITH EXTERNAL RECIRCULATION :


Take out the connection tube, and connect the in and out liquid tube as in the photo.

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Overtemperature alarm
Undertemperature alarm
Minimum working temperature
Maximum working temperature (depending on model).
Scale selection^º
External probe^º
Calibration

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STARTINGUP:

1. Fill the tank with liquid.
2. Connect the thermostat into the mains.
3. Turn on the main switch (2). For F(bath switch on the cooling unit if you wish to work with temperatures below the ambient. For other baths switch on the external cooling unit. (If necessary).

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4. Select the temperature (SV) with (5) and (6) buttons.
5. Press (7) to start up. While working, PV:(Point value) Shows the real fluid temperature and SV:(Set point) Shows the desired temperature.

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ACCESS T O CONFIGURATION PARAMETERS:

OVERTEMPERATURE

- Press (7) to stop the thermostat (STATUS: STOP).
- Press (3) one time, on the display appears: **"OVERTEMP"**, by pressing (5) or (6) the overtemperature degrees are selected. The alarm is activated when the read temperature (PV) exceeds the overtemperature degrees of the selected temperature (SV). The heating is cut off.
- The following screen appears:

where PV is the present temperature read.

- To remove the screen indications press (7), let cold the fluid and press again (7) to start.

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UNDERTEMPERATURE

- Press (7) to stop the thermostat (STATUS: STOP).
- Press (3) one time, on the display appears: **"UNDERTEMP"**, by pressing (5) or (6) the undertemperature degrees are selected. The alarm is activated when the read temperature (PV) is below the undertemperature degrees of the selected temperature (SV). The heating is cut off.
- The following screen appears:

where PV is the present temperature read.

- To remove the screen indications press (7), let cold the fluid and press again (7) to start.

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MINIMUM WORKING TEMPERATURE

- *Press (4) until the display shows **LOWER LIM**.*
- *By pressing (5) and (6) select the desired value.*
- *Press (3) until the display shows the initial screen.*

MAXIMUM WORKING TEMPERATURE

- *Press (4) until the display shows **LOWER LIM**.*
- *Press (3) until the display shows **UPPER LIM**.*
- *By pressing (5) and (6) select the desired value.*
- *Press (3) until the display shows the initial screen.*

WORK WITH FARENHEIT OR CELSIUS SCALE


- *Press (4) until the display shows **LOWER LIM**.*
- *Press (3) until the display shows **SELECT**.*
- *By pressing (5) and (6) select the desired scale.*
- *Press (3) until the display shows the initial screen.*

EXTERNAL PROBE

- *Press (4) until the display shows **LOWER LIM**.*
- *Press (3) until the display shows **EXT SENSOR**.*
- *By pressing (5) and (6) select YES or NO.*
- *Press (3) until the display shows the initial screen.*

CALIBRATION

- *Press (4) until the display shows **LOWER LIM**.*
- *Press (3) until the display shows **CALIBRATE**.*
- *By pressing (5) and (6) select the desired value.*
- *Press (3) until the display shows the initial screen.*

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ALARMS THAT CAN APPEAR ON THE DISPLAY:**OVERTEMPERATURE: OVERTEMP** (See page 7).**UNDERTEMPERATURE: UNDERTEMP** (See pag 7).**LEVEL:**

The screen shows:

LEVEL	STATUS
FAIL	ALARM

- 1- Cause: Low fluid level. Effect: Cut heating off. Action to solve:
1) Stop pressing (7). 2) Add fluid.
3) Start pressing (7).
- 2- Cause: The mechanical thermostat starts
Effect: Cut heating off.
Action to solve:
1) Stop pressing (7).
2) Reset the safety thermostat placed at the rear part of the machine. Review the temperature adjustment of the thermostat.
3) Start pressing (7).

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Effect: Cut heating off.

Action to solve:

- 1) Stop pressing (7).
- 2) Disconnect from the mains.
- 3) Take it to the Technical Service.

2- The screen shows:

RTD FAIL	STATUS
INT SHORT	ALARM

Cause: Shorted external temperature probe.

Effect: Cut heating off.

Action to solve:

- 1) Stop pressing (7).
- 2) Disconnect from the mains.
- 3) Take it to the Technical Service.

3- The screen shows:
FAIL

RTD	STATUS
EXT OPEN	ALARM

Cause: External temperature probe cut or disconnected.

Effect: Cut heating off.

Action to solve:

- 1) Stop pressing (7).
- 2) Check that the external probe is correctly connected.
- 3) Start pressing (7).

If the failure persists: 1)

- Stop pressing (7).
- 2) Disconnect from the mains.
- 3) Take it to the Technical Service.

4-The screen shows:

RTD FAIL	STATUS
EXT SHORT	ALARM

Cause: Shorted external temperature probe.

Effect: Cut heating off.

Action to solve:


- 1) Stop pressing (7).
- 2) Check that the external probe is correctly connected.
- 3) Start pressing (7).

If the failure persists: 1)

- Stop pressing (7).
- 2) Disconnect from the mains.
- 3) Take it to the Technical Service..

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External Probe connection:

The apparatus can be used to cool external equipment through a liquid circulation and an external probe. To install the water circulation circuit see section "INSTALLATION". The connector of the external probe is placed at the rear part of the apparatus head and the contacts are disposed according to figure "1". To read through this probe, select "EX SENSOR" according to the following configuration diagram:

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GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use of causes beyond the control of FUNGILAB, S.A.

Any manipulation of the apparatus by unauthorized personnel by FUNGILAB, S.A., cancels the guarantee automatically.

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- 9600 bauds, 7 dat bits, 1 bit stop, without parity no handshake.
- Up to 15m of cable.
- Cable connection DTE to DTE 9 pins male.
- Protocol ASCII with LRC against transmission errors
- MASTER-SLAVE Communication (PC computer as a master).
- The device never sends nothing if computer does not ask for it.

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- *Messages starts with character #.*
- *Device address must be included. This field goes from 00 to 99, size always 2 bytes in decimal-ASCII format.*
- *Data field is the message*
- *LRC is a calculated value to check good transmission.*
- *At the end of the frame always a CR (0D) followed by LF (0A).*
- *Slave device answers always sending its address*

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Frame includes a transmission error checking field called LRC. This field size 2 bytes. This field is calculated by the sender device and includes it on the frame. The receiver device also calculates and if are equal the transmission was good.

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LRC calculation procedure: 1.

1. Convert each character on it hexadecimal value.
2. Add
3. Take the smallest byte (two first characters starting from right hand side)
4. Subtract from FF on previous result
5. Add 1

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It could be found is any specilized purchaser:
Cable Null - modem de 9 pins

PC	Device
1	7, 8
2	3
3	2
4	6
5	5
6	4
7	1
8	

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Connectors in both cable sides are 9 pins female

INSTRUCCIONES RS- 232 EQUIPO FRIGITERM - ULTRATERM / RS-232 INSTRUCTIONS FOR FRIGITERM - ULTRATERM

PVT? Solicitar temperatura de la sonda PV (Point Value) *Get probe temperature PV (Point value):*
 Comando / *Command:* #01PVT?KK<CR><LF>
 Respuesta / *Answer:* #01+000,0KK<CR><LF>

SVT? Solicitar temperatura de consigna SP (Set point) *Get Set point temperature:*
 Comando / *Command:* #01SVT?KK<CR><LF>
 Respuesta / *Answer:* #01+000,0KK<CR><LF>

SVT Establece temperatura de consigna SP (Set point) *Set point temperature*
 Comando / *Command:* #01SVT<BL>+000,0KK<CR><LF>
 Respuesta / *Answer:* #01OKKK<CR><LF>

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MAINTENANCE

Before removing the side covers disconnect the apparatus from the mains.

If tap water is used, there will be calcareous depositions on the bottom of the tank and on the heating element. Use of the heating element must be avoided when covered with these depositions because it loses calorific power and it is destroyed quicker.

It is advisable to clean the inner tank and heating element periodically with a proportion of 10% of hydrochloric acid or with a solution of water and vinager (acetic acid). Then rinse it with clean water.

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