## Noria R

## Rotary mixer adjustable



Rotary mixer adjustable, controlled by microprocessor.

- · For rotary motion with various · Allows place and removal of tubes from 6 to 35 mm diameter.
- · Adaptable to any tube dimension depending on the needs.
- Rotation speed adjustable, controlled by microprocessor.
- the tubes without stops the machine.
- Completely adjustable slope from 0° to 90°.
- Smooth and progressive start and stop.
- · Security clutch allows to stop rotary motion with the hand without problems.
- · Anti-scratching surface with chemical resistant coating.
- Suitable for continuous running.
- Packaging: 47x37x32cm



## **Technical Specifications:**

| Model             | NR50E       |  |
|-------------------|-------------|--|
| Reference         | 10000-01062 |  |
| Mixing            |             |  |
| Power (W)         | 50          |  |
| Type of movement  | Rotary      |  |
| Angle of slope    | 0°-90°      |  |
| Speed Range (rpm) | 5-50        |  |
| General Data      |             |  |
| Dimensions (mm)   | 320x260x200 |  |
| Weight (kg)       | 4,7         |  |
| Protection IP     | IP54        |  |

| Accessories                          | Reference   |
|--------------------------------------|-------------|
| Platform with clamps for:            |             |
| 25 tubes Ø 6-7 mm                    | 20000-00093 |
| 25 tubes Ø 8-9 mm                    | 20000-00094 |
| 25 tubes Ø 10-11 mm Eppendorf        | 20000-00095 |
| 20 tubes Ø 12-13 mm                  | 20000-00133 |
| 20 tubes Ø 14-15 mm                  | 20000-00134 |
| 20 tubes Ø 16-17 mm                  | 20000-00135 |
| 19 tubes Ø 18-19 mm                  | 20000-00136 |
| 16 tubes Ø 20-21 mm                  | 20000-00137 |
| 14 tubes Ø 22-23 mm                  | 20000-00138 |
| 14 tubes Ø 24-25 mm                  | 20000-00139 |
| 14 tubes Ø 27-30 mm                  | 20000-00140 |
| 10 tubes Ø 32-35 mm                  | 20000-00141 |
| Racks for:                           |             |
| 60 culture tubes Ø 12-13,5 mm        | 20000-00051 |
| 36 centrifuge tubes 15ml Ø 16-17,5mm | 20000-00052 |
| 20 centrifuge tubes 50ml Ø 29-31mm   | 20000-00053 |
| 60 tubes Eppendorf de 0,5 a 2 ml     | 20000-00096 |



